


4.0mm Diameter						5.0mm Diameter						6.5mm Diameter				7.5mm Diameter	
0°	15°	25°	0°	15°		0°	15°	0°	15°	25°	0°	15°	0°	15°	0°	15°	


SHORT IMPLANTS

3.0mm Post




MAX 2.5 MAXILLARY ANTERIOR IMPLANTS

2.5mm Post




NARROW IMPLANTS

2.0mm Post



RESTORATIVE/LABORATORY KIT



NOTES

- Snap-on sleeves are only specific for abutment diameter.
- Abutment height is not a criterion for proper selection of snap-on sleeves.
- Transfer dies correspond to exact diameter **and** height of abutment placed.
- Because of machining tolerances, sleeves may not reach the height of contour for some angled abutments.